


How to Connect to a lab in NetSupport

1. Open NetSupport




2. Select New




The screenshot shows the 'New Class' dialog box in NetSupport. The 'Lesson Details' section includes fields for 'Teacher username/email' (set to 'akruger'), 'Teacher name', and 'Lesson title'. There are also dropdown menus for 'Objectives' and 'Outcome', and a 'What time does this Lesson finish?' section with a clock icon and dropdowns for '9' and '15'. A checkbox for 'Open lesson' is checked. The 'Connect to' section displays a grid of icons for various labs: '302', '215', 'Lab 1', 'lab 2', 'Colby', and 'Eval'. Below these are icons for 'test' and 'TEST - IT'. A vertical toolbar on the right contains buttons for 'Select', 'New', 'Edit', 'Remove', and 'Power On'. The 'New' button is circled in red. At the bottom, there are checkboxes for 'Create a Journal' and 'Create a Student Register', and buttons for 'Configure', 'Quit', 'OK', and 'Help'.

3. Select Room Mode



4. Enter the name of the lab. Examples: 215, 302, LMC, MobileCart 1, MobileCart 2





5. Select your newly entered room then hit the Select button then hit OK

New Class ✕

Lesson Details

Teacher username/email: Teacher name: Lesson title:

Objectives: Outcome: What time does this Lesson finish?

Open lesson 2

Connect to

302 1 215 Lab 1 lab 2 Colby Eval

LMC test TEST - IT

Create a Journal Create a Student Register

3

Select
New
Edit
Remove
Power On

OK